SERVICE INTRODUCTION

Air volume — Air velocity — Particle contamination

Radel&Hahn Clean Room Servicing do servicing and repairwork of laminar boxes since 1980. Our company widened the range of his services available from 1988: we do ventilation measurement, adjustment, maintenance, repair work of ventilation systems of clean-rooms.

During maintenance the tests of **laminar boxes** are performed in harmony with regulations of the standard **MSZ EN ISO 14644**.

Air volume measurement

This test aims to measure the average airflow velocity and airflow evenness and inlet air volume inside the clean space because through it air exchange rate (rate in every hour) can be calculated as well.

During the test we use TESTO airflow measurements and TSI flow hood balometer, which are calibrated in accredited laboratories

Airflow velocity test

During the test we measure the average airflow velocity and airflow consistency according to the conditions defined in the standards and according to technical description.

Measurement of particle contamination:

This test aims to define the air cleanliness in the equipment and to do classification analysis.

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Leakage test of Hepa and Ulpa filter system

The aim of the test is to certify, that the filters of high filtration efficiency were properly integrated into the devices and they are leakage-free furthermore filters are free from damages (small holes, deterioration of filter material or of the frame sealant) or prove that they are free from any leakage (leakage at filter frame sealant or hole on the frame). We let DEHS - aerosol in front of the filter and we check immediately its quantity at the other side of the filter and the frame.

We use DOP SOLUTIONS LTD. device families calibrated by Hungarian distributor

Saturation analysys of Hepa and Ulpa filter

This analysys of saturation and the physical condition of HEPA filter based on the measured pressure difference between the two sides of the filter. This analysys informs us about the filter clogging, which keeps on increase with the time. Filter clogging and the implied pressure difference between the two sides of the filter cannot reach the upper limit given by the manufacturer, because it can cause that the material of the filter disrupts.

During the test ΔP measuring instruments are used, which are calibrated in accredited laboratories.

Meas-

urement of particle contamination

During the test we measure particulate contamination of the air. Based on the result we get, we diagnosticate the air cleannes rate inside the cleanroom and we classify the cleanroom into certain cleanliness class.

We use LASAIR III 350L air particle counters and these measuring instruments are calibrated in Hungarian accredited laboratories

Monitoring of air flow direction:

The goal of this test is that we monitor the direction of the laminar air flow in cleanzone. The smoke produced by a smoke generator makes the air-flow visible during measurement and we make a video record about it. We hand over this record to our customer on an electronic device.

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Customer: Radó & Hahn Zrt.

Equipment Model: DOP SP200 DMS

Equipment Type: Aerosol Photometer

Equipment Manufacturer: DOP Solutions, UK

Serial Number: 4046863

Date of Calibration: 22 Nov 2017

Validity of Calibration: 12 months from Date of Calibration (hours if Blank there is no flow meter) if the hours presented the calibration period, from the hours between the validity.

TEST EQUIPMENTS	Instrument	Type	SN	Cal Date	Cal. Valid Date
Digital Multimeter	Vision 10200		1013532107	Sept. 2017	Sept. 2018
DOP Photometer	DOP touch		413154	Nov. 2017	Nov. 2017
National Mine Safety System	DOP-350		413157	Nov.	N/A

Test Conditions:

- The instrument has been tested at 20°C ±0.1°C, RH: NC, Atmospheric pressure
- 20 minutes of stillness was passed
- The calibration test has been carried out after a full service to the instrument
- The instrument is tested against a calibrated test cell. The test cell is tested and calibrated annually
- The instrument is tested using our in-house documented test procedure.
- The instrument is tested in filter against a Sierra Mass Flow Meter flow meter, calibrated annually.
- The uncertainty of the calibration is based on the uncertainty of the contributors.
- The calibrations were measured on a Matestech MTEB digital multimeter, which is tested and calibrated annually.
- The calibration procedure meets the requirements for an instrument to be used for testing by external customers, as determined using the adjustable aerosol concentration (FSL)
- Calibration was performed at PCT Ltd. Calibration Lab.

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Calibrated by: [Signature]

Signature: [Signature]

Pressure condition — Temperature — Humidity

Pressure condition test

This procedure serves for measuring air pressure between areas of cleanroom. This measurement aims to prove that air supply system is suitable for maintaining pressure level between premises of the cleanroom and between the cleanroom and the outer air.

Based on customer request we do air pressure difference test and the audit log

Temperature and humidity distribution analysis

The aim of the analysis to prove that ventilation system of the cleanroom is suitable for keeping the air temperature and air humidity within certain limit values (relative air humidity or expressed as dew point) which is specified by the customer at this definite place

Reset test

We prove during this test that cleanroom can restore its condition to meet the requirements of a defined purity class again, after it was exposed to particle contamination for a short period of time. Performing this test is not recommended in workstations, where air streams laminarily.

Measurement of recovery time of material sluice, initial flooding level

Measurement of recovery time of material sluice, clarification level

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Customer recommended calibration PCT recommended calibration

Cal Due date: Nov. 2018
 customer asset #: n/a

PCT certifies that the instrument listed above meets or exceeds all published specifications and has been calibrated using equipment and standards traceable to the USA National Institute of Standards and Technology (NIST). The procedure used to calibrate this instrument is documented in the current PMS instrument service manual unless otherwise stated.

PCT procedure used: Aerosol Particulate Counter ASTM F-328 ASTM F-468 ASTM F-469 ASTM F-649

— resolution of the calibration method 10.0 µm for ASTM F-468 rule

calibration performed at PCT service lab calibration performed at customer site

INSTRUMENT CONDITION AS LEFT:

in tolerance operational failure
 out of tolerance physical damage

out of tolerance description:

CALIBRATION STANDARD USED:

PARTICLE SIZE	ACCURACY	LOT#	QNTY	EXP. DATE	PARTICLE SIZE	ACCURACY	LOT#	QNTY	EXP. DATE
0.29µm	+16mm	18221	Dec2009	Sep2018	5.02µm	+4.0mm	18218	Jan2014	Aug2018
0.50µm	+4.0mm	18214	Jan2014	Aug2018	0.29µm	+4.0.04µm	18210	Mar2011	Oct2018

CALIBRATION EQUIPMENT, INSTRUMENT, SIN, CAL DATE, CAL QTY, DATE, AKKRED #

DISP. METER	INSTRUMENT	SIN	CAL DATE	CAL QTY	DATE	AKKRED #
DISP. METER	Metrohm MY30	01027706	Jan. 2017	Jan.	2018	
AIR FLOW METER	Serra Flow Meter	102831	Nov. 2016	Nov.	2017	
REF. INSTRUMENT	Laser II-110 LRI	99103	May. 2017	May.	2018	

SUPPORT EQUIPMENT, INSTRUMENT, SIN

REHUMIDIFIER SYSTEM	INSTRUMENT	SIN
REHUMIDIFIER SYSTEM	PRECIFR10	12396
FILTER	PARTISLUK 2 FFP	2385
LASER-POWERMETER	LASER-MATE-C	2385

CERTIFIED BY: Ivan Nagy CALIBRATION DATE: 17. Nov. 2017

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Customer recommended calibration PCT recommended calibration

Cal Due date: Dec. 2018
 customer asset #: n/a

PCT certifies that the instrument listed above meets or exceeds all published specifications and has been calibrated using equipment and standards traceable to the USA National Institute of Standards and Technology (NIST). The procedure used to calibrate this instrument is documented in the current PMS instrument service manual unless otherwise stated.

PCT procedure used: Aerosol Particulate Counter ASTM F-328 ASTM F-468 ASTM F-469 ASTM F-649

— resolution of the calibration method 10.0 µm for ASTM F-468 rule

calibration performed at PCT service lab calibration performed at customer site

INSTRUMENT CONDITION AS LEFT:

in tolerance operational failure
 out of tolerance physical damage

out of tolerance description:

CALIBRATION STANDARD USED:

PARTICLE SIZE	ACCURACY	LOT#	QNTY	EXP. DATE	PARTICLE SIZE	ACCURACY	LOT#	QNTY	EXP. DATE
0.29µm	+16mm	18221	Dec2009	Sep2018	5.02µm	+4.0.04µm	18210	Mar2011	Oct2018
0.50µm	+4.0mm	18214	Jan2014	Aug2018	0.29µm	+4.0.04µm	18210	Mar2011	Oct2018

CALIBRATION EQUIPMENT, INSTRUMENT, SIN, CAL DATE, CAL QTY, DATE, AKKRED #

DISP. METER	INSTRUMENT	SIN	CAL DATE	CAL QTY	DATE	AKKRED #
DISP. METER	Metrohm MY30	01027706	Jan. 2017	Jan.	2018	
AIR FLOW METER	Serra Flow Meter	102831	Nov. 2016	Nov.	2017	
REF. INSTRUMENT	Laser II-110 LRI	99103	May. 2017	May.	2018	

SUPPORT EQUIPMENT, INSTRUMENT, SIN

REHUMIDIFIER SYSTEM	INSTRUMENT	SIN
REHUMIDIFIER SYSTEM	PRECIFR10	12396
FILTER	PARTISLUK 2 FFP	2385
LASER-POWERMETER	LASER-MATE-C	2385

CERTIFIED BY: Ivan Nagy CALIBRATION DATE: 01. Dec. 2017

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**You may download the full referenclist from our web-site:
www.radel-hahn.hu**

radel & hahn zrt

H-4028 Debrecen, Kassai út 92.

Telefon: 0036 52448441 Fax: 0036 52415258

E-mail: info@radel-hahn.hu Web-site: www.radel-hahn.hu